

## PRODUCT DATA SHEET



### Product description

|                     |                                  |
|---------------------|----------------------------------|
| Catalogue Number    | 1-202-C4-2                       |
| Product description | AHN myTip® 1000 µL, blue, racked |
| Sales Unit          | case / 4 x 8 x 96 tips           |

### Product characteristics

|                                |                                                           |
|--------------------------------|-----------------------------------------------------------|
| Graduations                    | yes                                                       |
| Filter                         | no                                                        |
| Filter Specifications          | N/A                                                       |
| Low Retention                  | no                                                        |
| Sterilised                     | no                                                        |
| Autoclavable (max. °C, period) | yes, 121 °C, 20 min                                       |
| Purity                         | Free of RNase, DNase, Human DNA, Pyrogens, PCR Inhibitors |

### Size

|                   |          |
|-------------------|----------|
| Max. volume       | 1000 µL  |
| Length of product | 84.30 mm |

### Material/colours

|                   |                    |
|-------------------|--------------------|
| Product material  | polypropylene (PP) |
| Colour of product | blue               |

### Packaging

|                                   |                    |
|-----------------------------------|--------------------|
| Type of Packaging                 | Racked             |
| Pcs./Rackbox                      | 96 tips            |
| Pcs./Pack                         | 768 tips           |
| Pcs./Case                         | 3072 tips          |
| Dimension Case Carton (H x W x D) | 399 x 273 x 400 mm |

### Storage and Handling

|             |                       |
|-------------|-----------------------|
| Temperature | below 35°C            |
| UV light    | protect from UV light |

**Note:** Improper storage can initiate degradation, which results in colour changes and can have negative effects on the physical properties of the products.

**Safety:** The product is not classified as dangerous.

**Recycling:** The product is recyclable.

This document was sent online and is valid without signature

#### CEOs

Alexander Spector  
Chirag Shah  
Henrique Hosken  
Tejas Shah  
Torsten Buchwald

#### Company Address

AHN Biotechnologie GmbH  
Uthleber Weg 14  
D-99734 Nordhausen

#### Contact

T: +49 (0) 3631 65242-0  
F: +49 (0) 3631 65242-90  
E: info@cappahn.com  
W: www.cappahn.com

#### Bank Account

Commerzbank AG €/€  
Konto-Nr.: 6065155  
BLZ: 82040000  
BIC: COBADEFF820  
IBAN: DE16820400000606515500

#### Commercial Register

HRB 404781  
Amtsgericht Jena  
USt-IdNr. DE812717024

